INSTRUCTION FOR USE

APERTO OTW PACLITAXEL RELEASING HEMODIALYSIS SHUNT BALLOON DILATATION CATHETER

Approved
25.08.2022

DEFINITIONS

PTA: Percutaneous Transluminal Angioplasty

OTW: Over-The-Wire design

Proximal: The catheter part nearest to the operator

Distal: The catheter part nearest to the patient

Stenosis: Abnormal narrowing of vascular pathways

DESCRIPTION

The APERTO OTW is a balloon dilatation catheter with an Over-The-Wire (OTW) design and a semi-compliant inflatable balloon mounted on the distal part of the catheter near the tip. It has a biluminal shaft design. The one lumen (guidewire lumen) permits the use of a guidewire with a maximum outer diameter of 0.035 inch. The other lumen (balloon inflation lumen) is used for balloon inflation and deflation with sterile saline solution or a mixture of saline and contrast media (1:1). The balloon is coated with 3.0µg of Paclitaxel per mm² balloon surface. The antiproliferative drug is released during balloon inflation to the vessel wall in order to prevent or reduce recurrence of restenosis.

The proximal part of the catheter includes one female luer-lock port connected to the inflation lumen, and one female luer-lock for the guidewire lumen.

Recommended introducer sheath compatibility:

Sheath size	Balloon Ø	Balloon length
6 F	4.0 – 6.0 mm	20 – 60 mm
7 F	7.0 – 10.0 mm	20 – 60 mm
8 F	12.0 mm	20 – 40 mm

Catheter Length: 40, 60 or 80cm

The hydrophilic coating which covers the distal part of the APERTO OTW reduces friction, facilitating passage through the vessels (see figure 1).

Figure 1: The hydrophilic coated part of the APERTO OTW

Two radiopaque markers placed in the cylindrical part of the balloon (one placed in the proximal part of the balloon and another one placed in the distal part) aid in positioning the balloon in the desired location. They, in conjunction with fluoroscopy, enable accurate positioning of the balloon within the vessel.

The APERTO OTW has a tapered and atraumatic tip for easy entry into peripheral arteries and to improve crossing of tight stenosis.

The APERTO OTW is available in different balloon sizes. Nominal balloon diameter and length are printed on the hub, on the pouch label and on the cardboard box label.

INDICATIONS FOR USE

APERTO OTW is indicated for Percutaneous Transluminal Angioplasty (PTA) in the peripheral vasculature, especially dialysis shunts and for treating obstructive lesions of native or synthetic arteriovenous dialysis fistulae. APERTO OTW can be also used for PTA of stenosis at subclavian vein, veins of the arm and the head, forearm vein, brachio-cephalic vein, vena cava superior and inferior.

CONTRAINDICATIONS

- Inability to cross-target lesion with a guide wire.

- Coronary arteries, arteriae pulmonales, aorta ascendens, arcus aorta, aorta descendens to the bifurcatio aortae, arteriae coronariae, arteria carotis communis, arteria carotis externa, arteria carotis interna, arteriae cerebrales, truncus brachiocephalicus, venae cordis, venae pulmonales, or the neuro-vasculature.
- Allergy, intolerance or hypersensitivity against Paclitaxel.
- Patients who exhibit angiographic evidence of existing thrombus
- Patients with contraindication for anti-platelet/anticoagulant therapy, including allergy.
- Known allergy to contrast media
- Pregnant women or women who are intending to become pregnant or men wishing to preserve their fertility
- Resistant (fibrotic or calcific) lesions, which cannot be predilated
- Total occlusion of the target vessel

POTENTIAL COMPLICATIONS

The following complications may occur as a result of PTA:

Puncture related:

- Local haematoma
- Local hemorrhages
- Thrombo-embolisation
- Thrombosis
- Arteriovenous fistula
- Spurious aneurysm
- Infections

Dilatation related:

- Dissection in the dilated vessel wall
- Perforation of the vessel wall
- Failure of prosthesis
- Prolonged spasms
- Acute re-occlusion with consequent surgical intervention
- Re-stenosis of the dilated vessel

As in every angiography, there may be an allergic reaction to the contrast medium.

PRECAUTIONS

Percutaneous transluminal angioplasty may only be carried out with this balloon catheter by interventionalists who have sufficient experience with transluminal interventional techniques. They should also have a thorough knowledge of current literature about treatment of the acute complications that may occur during PTA.

This balloon catheter is supplied in a two-part, sterile, tear-open package. Before the catheter is used, it is important to ensure that the packaging is intact. Catheters that have been in open or damaged packaging must not be used.

The general technical requirements for catheter insertion must be observed at all times. This includes, for example, flushing the components to be inserted in the vessel system with sterile, isotonic saline solution and the usual prophylactic, systemic heparinisation.

Caution: Replace the air in the balloon, the balloon lumen and the central lumen with liquid before the catheter is used.

Allergic reactions to contrast medium should be identified before PTA.

In order to activate the hydrophilic coating, it is recommended to wet the APERTO OTW with a sterile saline solution immediately prior to its insertion into the body. **CAUTION! Do not wet the balloon in order to protect the drug coating!**

Do not advance the catheter against significant resistance. The cause of this should be determined first e.g. via fluoroscopy.

Each catheter should be stored in a dark, cool and dry place.

WARNINGS

APERTO OTW is intended for single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Every catheter has to be inspected prior to procedure to verify the functionality and to make sure that neither the outer nor the inner package is damaged or opened. The catheter must not be used after the "Use By" date specified on the package.

Use only a mixture of contrast medium and saline solution (1:1) to fill the balloon. Do not use air or any gaseous medium to inflate the balloon!

Do not expose the catheter to organic solutions e.g. alcohol! Do not move the guidewire during inflation of the balloon!

The balloon must be completely deflated before removing it from the artery.

Do not use a second drug-coated balloon at the same treatment site. Additional APERTO OTW balloon catheter can be used to treat the lesion longer than the maximum balloon length available under the following circumstances:

- each individual segment should be treated only once with a drug-coated balloon
- try to minimize overlapping of treated segment
- treat each segment with a new PTA balloon

Overlapping of consecutively placed balloons with a treated segment should be avoided by precisely angiographically positioned catheter, using the marker bands.

In order to eliminate as far as possible the balloon rupturing during use, the **Rated Burst Pressure (RBP)** must **NEVER BE EXCEEDED**.

PREPARATION OF THE CATHETER

First of all, it is very important to check that the packaging is neither damaged nor opened before use.

The catheter and balloon in the deflated state contain small amounts of air that must be eliminated before use.

Replace the air in the balloon, the balloon lumen and the central lumen with liquid before the catheter is used. The following steps should be taken:

- Flush out the guidewire lumen with heparinised saline solution.
- Attach a syringe partially filled with a sterile mixture of saline solution and contrast medium to the balloon tube.
- Remove the balloon protective sleeve.
- Whilst the syringe and balloon are held down, draw the air out of the balloon into the syringe. Repeat until there is absolutely no air left in the balloon. Repeat this process several times if necessary.
- Release the plunger and allow the saline to fill the folded balloon. **DO NOT** inflate the balloon.
- Repeat the two last steps until all air is removed from the balloon catheter.
- Replace the syringe with a dilatation syringe and a manometer that have been prepared beforehand, i.e. partially filled with saline solution.
- Draw the saline solution out of the catheter into the dilatation syringe

CATHETER INSERTION

The APERTO OTW is introduced percutaneously via a surgical cut or through an appropriate introducer sheath. Use of an introducer sheath is generally recommended.

Under fluoroscopic control the balloon should be pushed through the appropriate sheath over the guidewire and positioned at the lesion to be dilated.

The recommended balloon inflation and dilatation time is 90 sec.

CATHETER REMOVAL

APERTO OTW is removed after PTA over the positioned guidewire and through the introducer sheath. As the balloon may have been dilated several times, there may be some resistance when it is removed.

STORAGE

- Store at controlled room temperature in a dry, cool and dark place.
- Use the catheter before to its "Use by Date" date printed on the package label.
- **Do not** use the catheter if the packaging is opened or damaged.
- **Do not** expose the catheters to organic solvents, ionising radiation or ultraviolet light.

STERILISATION

This catheter has been sterilized with ethylene oxide gas. It will remain sterile for as long as the packaging remains unopened and undamaged.

Warning: For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

WARRANTY

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the APERTO OTW as described in this publication. Under no circumstances shall CARDIONOVUM be liable for any specific law. No person has the authority to bind APERTO OTW to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

CARDIONOVUM will not be responsible for any direct, incidental, or consequential damages resulting from re-use of the product.


Explanation of symbols used on the package labels.

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Cardionovum GmbH, Am Bonner Bogen 2, 53227 Bonn, Tyskland