RESTORE PLUS DCB PACLITAXEL RELEASING PTCA BALLOON CATHETER

Description

RESTORE PLUS DCB PTCA catheter is a double Lumen coronary balloon dilatation catheter for rapid exchange, with a balloon that can be enlarged and two radiopaque markers (proximal and distal), to facilitate the balloon positioning under X-ray fluoroscopy. The RESTORE PLUS DCB PTCA balloon is coated with 3.0 µg Paclitaxel per mm². The dilatation balloon is designed to inflate to a known diameter and length at recommended inflation pressures (refer to Compliance Chart). The proximal end of the RESTORE PLUS DCB PTCA catheter has a single female Luer-Lock port to inflate/deflate the balloon. The RESTORE PLUS DCB PTCA is a Rapid Exchange catheter with a working length of 140 cm and a guide wire lumen of 25 cm length. The RESTORE PLUS DCB PTCA catheter is compatible with guide wires of 0.014" (0.36 mm) diameter and guiding catheters with 5F.

Packaging

Delivered in a peelable pouch and cardboard box.

One unit per box (one catheter with a protection sheet and a protection wire in a dispenser). Sterilized with Ethylene Oxide Gas.

Do not use if the package is opened or damaged. Use before the expiry date clearly indicated on the label.

Storage: Keep in a dry place between 0-25°C. Keep away from light.

Indications

The application of the RESTORE PLUS DCB PTCA catheter is indicated for patients with symptomatic ischemic heart disease because of discrete lesions of the coronary arteries with a reference vascular diameter of 1.5 mm to 4.00 mm, in order to improve the coronary lumen diameter. The RESTORE PLUS DCB PTCA catheter preferentially is applied for the dilation of coronary in-stent restenosis and for dilatation of small diameter coronary artery lesions, below 3.00 mm.

Contra-Indications

- Known sensitivity to Paclitaxel.
- Patients who are not candidates for coronary bypass surgery.
- Patients with totally obstructed coronary arteries.
- Patients who exhibit angiographic evidence of existing thrombus.
- Patients with contraindication for anti-platelet/anticoagulant therapy, including allergy. This includes patients who have had major surgery, an obstetrical delivery, organ biopsy, or puncture of a non-compressible vessel within 14 days of this procedure. Also excluded are those patients with a history of gastrointestinal bleeding, recent C.V.A., diabetic hemorrhagic retinopathy, of any other condition compromised by prolonged anti-coagulation.
- Pregnant women or women of childbearing potential.
- Target lesions distal to a 50 % or greater stenosis, which cannot be predilated or target lesions proximal to untreatable areas of significant flow compromising disease.
- Resistant (fibrotic or calcific) lesions, which cannot be predilated (lesions resistant to complete balloon inflation at 20 bar).

Warning

- Touching of the balloon or contact with liquids prior to insertion should be strictly avoided, as this could damage the balloon coating.
- A pre-dilatation with a conventional balloon in complex lesions is recommended. An excessive use of the RESTORE PLUS DCB can lead to the detachment of the coating.
- Thrombocyte aggregation inhibitor therapy (antiplatelet drug therapy) (75 mg Clopidogrel + 100 mg Aspirin over 2 Months).
- Connect only Luer lock fitting to the hub of the RESTORE PLUS DCB PTCA catheter.
- Should unusual resistance be met at any time during lesion access, first try to remove the balloon and then, if not working, try to remove both balloon catheter and guiding catheter as a single unit.
- Check that the sterile pouch is not damaged before use. Do not re-sterilize or re-use, destroy product after use.
- When the balloon catheter is exposed to the vasculature, it should be manipulated while under high-quality fluoroscopic observation.
- No very long term studies have been performed to evaluate the carcinogenic or genotoxic potential of the drug paclitaxel coated coronary balloon and there are no adequate and well-controlled studies published in women and men intending to have children or pregnant women. Therefore, we do not know at the present time whether paclitaxel is safe to be used in any kind of treatment during pregnancy or in patients intending to have children.

To reduce the potential for damage in coronary vessels and to avoid any balloon rupture, the balloon pressure should not exceed the rated burst pressure (RBP cf. compliance chart).

Precautions for use

- The cardiology team must be sufficiently familiar with the products and their reference systems to avoid errors in the selection of equipment.
- Only physicians who have received appropriate training should use the device.
- Catheter placement should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
- Ensure that the medical team is trained on products and their reference system to avoid any error in choosing equipment.

After use, eliminate the product according to safety requirements related to the products contaminated by blood.

Potential Adverse Effects / Complications

Complications that may result from this procedure include (Complications can directly result in the patient's death):

- Possible allergic reactions to Paclitaxel
- Acute myocardial infarction
- Arterial perforation or rupture
- Bleeding complications or haematoma
- Coronary spasm; coronary embolism; coronary thrombosis
- Dissection or perforation of the coronary artery
- Injury or rupture of the coronary artery
- Total occlusion
- Ischemia
- Unstable angina
- Disturbance of cardiac conductivity
- Restenosis of the dilated artery
- Ventricular fibrillation
- Arterial spasm

Vascular events:

- Perforation or dissection of the coronary artery or aorta
- Cardiac perforation
- Emergency cardiac surgery
- Aneurysm formation
- Iatrogenic arteriovenous fistula formation
- Femoral artery dissection at vessel access site
- Vessel recoil and restenosis
- Thrombosis or occlusion
- Peripheral ischemia
- Distal embolization (air, tissue debris, thrombus)
- Local inflammatory reactions
- Resulting from possible injuries of the vascular wall

Balloon catheter events:

- Failure to reach or cross lesion
- Inflation difficulties
- Balloon rupture or pinholes
- Balloon deflation difficulties

Instructions for use

Prior to using the RESTORE PLUS DCB PTCA catheter from the package, inspect for bends, kinks, and other damage. Do not use if any defects are noted.

Choice of the RESTORE PLUS DCB PTCA catheter:

The balloon diameter, when the balloon has been inflated to its nominal pressure should not be larger than the artery proximal and distal to the stenosis.

Catheter preparation

Prepare the inflation device according to the manufacturer's instructions.

- Do not let the balloon come into contact with fluids prior to insertion.
- Do not allow the balloon to come into contact with the skin.
- Fill a 10 or 20 cc syringe with sterile saline solution.
- Remove the protection sheet and wire from the balloon of the catheter and discard them.

- Prepare a diluted solution of contrast medium and sterile saline solution, in accordance with the contrast medium manufacturer's instructions.
- Fill a 10 or 20 cc luer lock syringe with 1 to 3 cc of the diluted 1:1 contrast medium.
- Connect a stopcock to the hub of the RESTORE PLUS DCB PTCA catheter.
- Connect the syringe to the stopcock, Position the syringe vertically with the plunger pointing upwards. Apply strong suction with the syringe to evacuate all air from the balloon.
- Maintain the suction for 15 to 20 seconds and make sure that no more bubbles are seen passing through the diluted contrast medium. Release the plunger carefully and, if necessary, repeat the procedure, in case of residual bubbles.
- Maintain the vacuum in the RESTORE PLUS DCB PTCA catheter and close the stopcock.
- Connect the inflation system, avoiding any bubble presence. Open the stopcock and inflate the balloon to its nominal pressure (see compliance table). Verify the integrity of the system. In case of defects, do not use the catheter. After inflation, apply negative pressure and close the stopcock.

Catheter manipulation

Attach the "Y" connector to the guiding catheter, which has been previously placed in the vasculature. The choice of the guidewire and the "Y" connector remains the responsibility of the physician within a compatibility limit of 0.014" for the guidewire.

Insert the guidewire (0.014" max.) into the guiding catheter and advance it through the tip.

To avoid any movement of the guidewire, tighten the knob of the "Y" connector so that it closes firmly around the guidewire.

Back load the guidewire into the distal tip of the catheter ensuring that it exits through the notch located approx. 25 cm proximal to the dilatation catheter tip.

Open the valve and advance the catheter to the distal end of the guiding catheter. Caution: Always advance the catheter fully deflated and always on the guidewire. Continue the procedure according to accepted coronary angioplasty technique. The radiopaque markers aid the positioning of the balloon in the stenosis.

Balloon dilatation and Paclitaxel release

To ensure an effective release of Paclitaxel to the lesion site of the artery wall, the balloon should be inflated at the lesion with an effective balloon dilatation pressure between the nominal balloon inflation pressure up to the rated burst pressure, for the duration of at least 30-45 seconds. Do not exceed the rated burst pressure!

Catheter exchange procedure

Loosen the "Y" connector.

Hold the guidewire and the "Y" connector in one hand, while grasping the catheter balloon shaft in the other hand. Maintain the guidewire's position in the coronary artery by holding the wire stationary and begin pulling the dilatation catheter out of the guiding catheter.

Liability

Cardionovum has endeavoured to ensure that the devices comply with all relevant standards and regulations currently in force and to make sure that the quality of the devices meets the requirements of the above mentioned standards and regulations for a period upon the indicated expiry date. The above statement does not apply where the devices are used for a purpose other than its intended purpose.

Advice

The product is for single use only and should not be cleaned, disinfected and re-sterilized. This form for condition ing would lead to unsterile Products and could cause infections to the patient and also risk of balloon rupture.

Explanation of symbols used on the package labels.

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Cardionovum GmbH, Am Bonner Bogen 2, 53227 Bonn, Tyskland